Calibra® Ceram

Adhäsiver Kompositement

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- Calibra Ceram elect[†] Adhesive, PrimeBond[†] Adhesive or PrimeBond[†] Adhesive or PrimeBond[†] Adhesive
- Light curing of applied adhesive may be accomplished before or after seating restoration with Calibra Ceram Cement.
- Following application and/or drying (solvent evaporation), proceed immediately to cementation, dual cure cleanup and Light Transmissible Visible light curing through the restoration.

[†]Zirconium-based restorations are generally not light-transmittable. Cure adhesive on the preparation prior to cementation.

Technique Tip: The adjacent teeth and/or the external surfaces of the restoration may be lubricated with a water soluble medium to ease clean up of excess cement.

3.1 Crown/bridge/inlay/onlay/post cementation
Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditionning/Dentin Pretreatment/Adhesive Application. Refer to manufacturer's Directions for Use for selected adhesive. Adhesives (plus activators, if applicable) compatible with both light-cured, CO initiated methacrylates and self-cured methacrylates are recommended. Apply adhesive (plus activator if applicable) to the internal bonding surface of the restoration if directed by adhesive manufacturer's directions for use.

3.1.1 Cementation technique

- Danger of injury due to excessive force**
- Apply slow and steady pressure on the syringe
- Do not use excessive force – dual barrel syringe rupture may result

CAUTION

- Dual Barreled syringe dispensed**
- Remove syringe cap. Dispense and discard a small amount of the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully wipe away excess so base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement following use.
- Install a mixing tip on the cartridge by lining up the w-shaped notch on the outside of the mixing tip with the w-shaped notch on the syringe flange. Turn colored mixing tip cap 90 degrees in a clockwise direction to lock in place on syringe.
- Gently depress syringe plungers to begin the flow of material. DO NOT USE EXCESSIVE FORCE. If force is encountered, remove syringe from operating field, remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tips as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.

3.1.2 Apply a uniform layer of cement on the entire inner surface of the restoration directly from the mixing tip. For inlays/onlays, it may be helpful to apply a thin layer of cement to the internal portions of the tooth preparation to avoid any porosity or voids. **Technique Tip:** The mixing tip may be bent slightly to allow direct intraxial access for placement of cement into preparations with internal anatomy. For endodontic post spaces, use a Lentulo Spiral or metal file to aid placement in the post space is recommended. At room temperature, Calibra Ceram Cement requires a minimum work time of 2 minutes when protected from ambient light. If loaded into restoration lined with adhesive or adhesive/activator mixture, room temperature work time (restoration seating time) will be less (see Precautions).

3.1.3 Immediately seat the restoration in the mouth with gradual pressure. Verify complete seating. A gentle rocking or vibratory motion may be helpful to insure optimal seating.

- Incomplete cementation**
- Intralux work time is reduced due to warmth and contact with adhesive
- Assure restoration is fully seated within 1 minute
- Placement of multiple units may require multiple mixes of cement

3.1.4 Protect restoration from contamination and movement until the final set of the cement (5 minutes from start of mix or in the case of light-transmissible restorations, completion of light curing).

- Inadequate retention due to insufficient stabilization**
- Stabilize restoration during cleanup
- Protect restoration from movement during setting
- Allow cement to self-cure without disturbing for 5 minutes

3.2 Cleaning marginal excess Calibra Ceram Cement
All shades of Calibra Ceram Cement are radio-opaque, with a 1mm radio-opacity equivalent to 2mm radio-opacity of aluminum. Aluminum has a radio-opacity equivalent to that of dentin. This 1mm of material having a radio-opacity equivalent to 1mm of aluminum has a radio-opacity equivalent to that of dentin.

3.2.1 Self-cure cleanup
The excess cement will reach the “gelled” state after approximately 1-2 minutes in the mouth, allowing easy removal. Excess cement will remain in the “gelled” state for approximately 1 minute. If exposed to directed operatory light, “gali” state may be reached sooner and remain “gelled” for a shorter period. Immediately after reaching the “gelled” state, floss very lightly, only in the direction of restoration seating, to remove excess cement. Complete excess cement removal using an explorer with a rubber tip, a scaler or an explorer. **NOTE:** Cement within the crown has not yet set. Do not move, touch, or disturb the crown during cleanup.
Technique Tip: Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

3.2.2 Optional dual cure cleanup
Due to the dual cure property of Calibra Ceram dual barrel dispensed cement, the operator has the option of utilizing a curing light to facilitate cleanup. Light curing to facilitate cleanup must be accomplished within the first minute following intraxial insertion of the restoration. Excess cement cleanup may begin immediately using a brief exposure with the curing light. Conventional powered quartz tungsten halogen or LED lights producing only one peak wavelength around 470nm are recommended. Use of narrow spectrum output light in high power, broad spectrum output halogen or LED lights may produce unexpected results. Light curing mixed cement after one minute, or continuously at one spot for more than 5 seconds, at any time, will cause cement to adhesively set, making cleanup difficult.

Immediately after seating is verified, briefly light-cure excess cement at the margins by constantly moving the curing light tip around the margin for no more than 5 seconds per surface. Excess cement will reach a “gelled” state after this brief (< 5 seconds).

- Residual excess cement – overcure during cleanup**
- Morowave (470nm) LED or conventional light sources are recommended
- High power, dual or broad spectrum lights may cause premature hardening of excess cement
- Evaluate curing/light-cement interaction in laboratory prior to clinical use

The excess cement will remain in the “gelled” state for approximately 45 seconds following light exposure. All excess cement must be removed before final self-cure set is achieved, as outlined above. **NOTE:** Cement within the crown has not yet gelled or set. Do not move, touch, or disturb the crown during cleanup. Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

3.3 Curing, finishing and dismissal
Non-Light-Transmissible restorations (Self-Cure)
For metallic, thick or heavily opacified ceramic composite, or restorations that otherwise impede the transmission of light, once cleanup is completed and restoration is stabilized, allow Calibra Ceram Cement to self-cure without disturbing for 5 minutes from start of mix. Protect restoration from contamination and movement during the setting time. Following the self-cure set, check and adjust occlusion and polish as necessary. Patient may then be dismissed.

- Inadequate retention due to insufficient curing**
- Stabilize restoration during cleanup
- Protect restoration from movement during setting
- Allow cement to self-cure without disturbing for 5 minutes

Light-Transmissible restorations (Light-Cure)
For most non-metallic, light-transmissible ceramic or composite restorations, Calibra Ceram Cement may be visible light cured. **Important Technique Note:** If following special Instructions for In-office Milled Ceramic Restorations utilizing uncured adhesive technique, Calibra Ceram Cement must be visible light cured through the restoration after seating. Once cleanup is completed and restoration is stabilized, light cure all areas of the restoration using a visible light, curing unit (output spectral peak including 470nm, broad spectrum output light recommended), with a minimum output of 550mW/cm² or 20 seconds from each direction – buccal, lingual and occlusal. Following the light-curing, check and adjust occlusion and polish as necessary. Patient may then be dismissed.

- Inadequate polymerization due to insufficient curing**
- Check compatibility of curing light
- Check curing cure cycle
- Check curing output before each procedure

Important Technique Tips:

- When simultaneously cementing multiple single units or bridgework, it is recommended to employ the light-cure cleanup on one or two adjacent units only, allowing other units’ excesses to self-cure, providing ample cleanup time.
- Clean excess cement from metal instruments immediately as set cement will adhere to the instrument.

3.4 Finishing and polishing

- Removal of resin flash is best accomplished with the Enhancer’s Finishing System’s tips, cups and discs. The Enhancer’s System will remove flash and finish restoration margins without removal of or trauma to the enamel.
- Polish final restoration using Enhance[®] PoGo[®] Polishing system or Prisma[®] Gloss[®] Polishing Paste and Prisma[®] Gloss[®] Extra Fine Polishing Pastes. (See complete Directions for Use supplied with polishing product chosen)

4. HYGIENE
Cross-contamination – Infection

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Syringe cannot be reprocessed. Dispose of contaminated syringe in accordance with local regulations.

4.1 Cleaning and disinfection

- For dual-barreled syringe, remove used mixing tip and discard appropriately. Replace original syringe cap prior to storing. The dual-barreled syringe may be cleaned by scrubbing with a disposable towel soaked with hot water and soap or detergent.
- To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid package contamination. The use of protective barriers is an additional precautionary measure against gross contamination but not against all contamination.
- Contaminated debrief instruments should be cleaned and disinfected with a water-based hospital-level disinfection solution according to national/local regulations.
- Repeated disinfection may damage label.
- NOTE:** Destruction of the label by vigorous wiping. Wipe syringe gently.

5. LOT NUMBER AND EXPIRATION DATE
1. Do not use after expiration date. ISO standard is date: “YYYY/MM”
2. The following numbers should be quoted in all correspondence:

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