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3.3 Cementation technique		
	Dual Cure – Non-Light Transmissible Restorations	Light Cure – Light Transmissible Restorations
Type of Restoration	Self-cure	All
	PFM, Zirconia, Alumina, opaque ceramics & composites	Translucent ceramics & composites
1. Fill and seat crown	Up to 2 min	Up to 2 min
	Protect restoration from movement during gel phase cleanup through final set	
2. Remove Excess		
Stage from seating:	1-2 min	Up to 5 sec per surface/ light cure
	Gel Duration: 1 min	45 sec
3. Stabilize		
	Constant Occlusal Pressure	Light cure margins 20-40 sec
4. Final Set	6 min from start of mix	Up to 5 sec per surface/ light cure
	6 min from start of mix	Upon completion of light curing each surface/ 10 sec

1. Buccal and Lingual 2. Buccal, Lingual and Occlusal

- Danger of injury due to excessive force**
 - Apply slow and steady pressure on the syringe
 - Do not use excessive force – syringe rupture may result

3.3.1 Dual Barreled syringe dispensed

1. Remove syringe cap. Dispense and discard a small amount of material from the dual-barreled syringe. Be sure material is flowing freely from both ports. Holding syringe vertically, carefully flow away excess so base and catalyst do not cross contaminate and cause obstruction of the ports. Save syringe cap for replacement, following use.

2. Install a mixing tip on the cartridge by lining up the v-shaped notch on the outside of the mixing tip with the v-shaped notch on the syringe flange. Turn colored mixing tip cap 90 degrees in a clockwise direction to lock in place on syringe.

3. Gently dispense syringe plungers to begin the flow of material. DO NOT USE EXCESSIVE FORCE. If force is encountered, remove syringe from operating field, remove and discard mixing tip. Check for obstruction and confirm material is flowing from both syringe barrels. Wipe barrels and install new mixing tip as outlined above. Dispense a small amount through the mixing tip onto a mixing pad and discard.

4. Without delay, using gentle pressure, apply a thin, uniform layer of cement to the entire internal surface of the restoration directly from the mixing tip. At room temperature, Calibra Universal Cement offers a minimum work time of 2 minutes. **Technique Tip:** The mixing tip may be bent slightly to allow direct intraoral access for placement of cement into preparations with internal anatomy. For endodontic, post spaces, use a Lentulo Spiral or metal file to aid placement in the post space is recommended.

5. Immediately seat the restoration in the mouth. Verify complete seating. A gentle rocker or vibratory motion may be helpful to insure optimal seating.

6. Protect restoration from contamination and movement until the final set of the cement (6 minutes from start of mix or in the case of light-transmissible restorations, completion of light curing).

3.4 Cleaning Marginal Excess Calibra Universal Cement

3.4.1 Self-cure cleanup

The excess cement will reach the “gelled” state after approximately 1-2 minutes in the mouth, allowing easy removal. Excess cement will remain in the “gelled” state for approximately 1 minute. If exposed to directed operator light, “gel” state may be reached sooner and remain “gelled” for a shorter period. Immediately after reaching the “gelled” state, floss interproximally, only in the direction of restoration seating, to remove excess cement. Complete excess cement removal using an instrument such as a rubber tip, a scalpel or explorer. **NOTE:** Cement may remain “gelled” for a shorter period. Immediately after reaching the crown during cleanup. **Technique Tip:** Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

3.4.2 Optional dual cure cleanup

Due to the dual cure property of Calibra Universal Cement, the operator has the option of utilizing a curing light to facilitate cleanup. Light curing to facilitate cleanup must be accomplished within the first minute of cement placement. Intraoral curing may be achieved by the use of a curing light. **NOTE:** Cement within the crown has not yet cured or set. Do not move, torque, or disturb the crown during cleanup. **Technique Tip:** Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

Residual excess cement – overcure during cleanup

- Monowave (470nm) LED or conventional halogen lights are recommended
- High power, dual or broad spectrum lights may cause premature hardening of excess cement
- Evaluate curing light/cement interaction in laboratory prior to clinical use.

The excess cement will remain in the “gelled” state for approximately 45 seconds following light exposure. Light-cure cleanup must be completed before final set is achieved, as outlined above. **NOTE:** Cement within the crown has not yet cured or set. Do not move, torque, or disturb the crown during cleanup. Following all excess removal, exposed margins may be light cured 20-40 seconds to assist restoration stabilization.

3.4.3 Inadequate retention due to insufficient stabilization

- Stabilize restoration during cleanup.
- Protect restoration from movement during setting
- Allow cement to self-cure without disturbing for 6 minutes

Light-Transmissible Restorations (Self-Cure)
For metallic, thick or heavily opaque ceramic or composite, or restorations that otherwise impede the transmission of light, once cleanup is completed and restoration is stabilized, allow Calibra Universal Cement to self-cure without disturbing for 6 minutes from start mix. Protect restoration from contamination and movement during the setting time. Following the self cure set, check and adjust occlusion and polish as necessary. Patient may then be dismissed.

3.5 Curing

Non-Light-Transmissible restorations (Self-Cure)
For metallic, thick or heavily opaque ceramic or composite, or restorations that otherwise impede the transmission of light, once cleanup is completed and restoration is stabilized, allow Calibra Universal Cement to self-cure without disturbing for 6 minutes from start mix. Protect restoration from contamination and movement during the setting time. Following the self cure set, check and adjust occlusion and polish as necessary. Patient may then be dismissed.

Inadequate retention due to insufficient stabilization

- Stabilize restoration during cleanup.
- Protect restoration from movement during setting
- Allow cement to self-cure without disturbing for 6 minutes

Light-Transmissible Restorations (Light-Cure)
For metallic, light-transmissible ceramic or composite restorations, Calibra Universal Cement may be visible light cured. Once cleanup is completed and restoration is stabilized, light cure all areas of the restoration using a visible light, curing unit, designed to cure CO initiated methacrylates (spectral output including 470nm), with a minimum output of 550mW/cm² for 10 seconds from each direction (buccal, lingual and occlusal). Following the light curing, check and adjust occlusion and polish as necessary.

- Inadequate polymerization due to insufficient curing**
 - Check compatibility of curing light.
 - Check curing cycle.
 - Check curing output before each procedure

3.6 Finishing and polishing

Removal of resin flash is best accomplished with the Enhance[®] Finishing System of points, cups and discs. The Enhance System will remove flash and finish restorations without removal of or trauma to the enamel. Polish final restoration using Enhance[®] PoGo[®] Polishing system or Prisma[®] Gloss[®] Polishing Paste and Prisma Gloss Extra Fine Polishing Pastes. (See complete Directions for Use supplied with polishing product chosen). All shades of Calibra Universal Cement are radio-opaque, with a 1 mm radio-opacity equivalent to 1.8 mm radio-opacity of aluminum. Aluminum has a radio-opacity equivalent to that of dentin. Thus 1 mm of material having a radio-opacity equivalent to 1 mm of aluminum has a radio-opacity equivalent to that of dentin.

4. HYGIENE

Cross-contamination Infection
Do not reuse single use products. Dispose of in accordance with local regulations. To prevent the syringes from spilling or spurt of body fluids, use the following guidelines:

- Hands it is mandatory that the syringes are handled with clean/disinfected gloves. Do not reuse syringes if contaminated.
- Syringe cannot be reprocessed. Dispose of contaminated syringe in accordance with local regulations.

To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended. The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination.

Incidental contact of the syringe with water, soap or a water-based hospital-level disinfection solution will not damage syringe body. Do not allow any solution contact with contained material. Discard composite material that has been in contact with any fluid or non-sterile instrument.

NOTE: Destruction of the label by vigorous wiping. Wipe syringe gently.

5. LOT NUMBER AND EXPIRATION DATE

1. Do not use after expiration date. See standard uses: “YYMM-MM-DD.” The following numbers should be quoted in all correspondences:

- Reorder number
- Lot number
- Expiration date

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