



CEREC® CAD/CAM Solutions

Designed to simply work better together

Celtra® Duo (ZLS) blocks, Prime&Bond Elect® adhesive, and Calibra® Ceram cement were designed to enhance and strengthen the individual benefits each of them provides, resulting in an easy-to-use system that streamlines the restoration process.

Celtra Duo (ZLS) blocks

- Restoration longevity of Celtra Duo (ZLS) is ensured when used with Prime&Bond Elect adhesive and Calibra Ceram cement
- Firing is optional: choose either fire and seat or polish and seat

Prime&Bond Elect Adhesive

- No need to use a self cure activator when used with Calibra Ceram cement
- Low film thickness to allow passive seating of the crown

Calibra Ceram Cement

- One-step curing when used with Prime&Bond Elect adhesive
- 10-second tack cure window and 45-second gel phase ensures an easy, no-stress cleanup



The Celtra® Duo (ZLS) 3-Step Restoration and Cementation System for CEREC® Users Designed to simply work better together

STEP 1

Design and mill the restoration

Design the restoration as usual with CEREC, then mill it out using Celtra Duo (ZLS) material.

Celtra Duo (ZLS) advantages:

- You choose the processing pathway: fire and seat, or polish and seat—you're always in control



STEP 2

Ⓐ Apply Prime&Bond Elect® Adhesive

Apply Prime&Bond Elect adhesive to the tooth.

Prime&Bond Elect adhesive advantages:

- Universal application means you're in control: self-etch, total-etch and selective-etch—it's always your choice
- Low film thickness
- Virtually no post-op sensitivity



Ⓑ Apply Calibra® Ceram Cement

After etching and silanating the intaglio surface of the restoration, apply a thin, uniform layer of Calibra Ceram cement to the internal surface of the restoration.

Calibra Ceram cement advantages:

- High bond strength for long-term restoration success
- Easy excess cement cleanup:
 - wide tack cure window of up to 10 seconds means no worry of over-curing
 - 45-second extended gel phase gives you the time you need for a thorough and effective cleanup



STEP 3

Seat restoration

Seat the restoration; the cement will set permanently after final light curing of all areas of the restoration.

Dentsply Sirona Restorative
38 West Clarke Avenue / Milford DE 19963
800-532-2855 / www.dentsply.com

©2016 Dentsply Sirona Inc. All rights reserved.

MLA-51N000 (8/3/16) 

Developed to Make a Difference.

To place an order, contact your Dentsply Sirona Representative or call 855-723-5872.