- Attach supplied needle tip to end of the Calibra® Silane Coupler syringe. Gently apply pressure to syringe plunger. DO NOT USE EXCESSIVE FORCE. Silane agent should express easily, one drop at a time. Directly apply the silane to the clean internal surface of the restoration. Air-dry the silanated surface of the restoration with clean, dry air from a dental syringe or allow to air dry. **NOTE:** Avoid application of the silane to the external surface of the restoration. Contact with the external surface will cause bonding of the cement to the external surface, complicating clean-up and necessitating removal of the glazed portion of the external surface of the restoration.
- If silanated surface becomes contaminated, clean surface with alcohol or acetone, allow to air dry and repeat application as outlined above.

3.3 Tooth conditioning/dentin pretreatment/adhesive application

- Following restoration try-in, rinse preparation thoroughly with water spray and air dry.
- Proceed with adhesive application as directed in the complete directions for use.

<div>Delayed or reduced adhesion due to adhesive incompatibility</div> <div> <ul style="list-style-type: none">Select compatible light cure adhesive Follow selected adhesive directions for use</div>	<div>CAUTION</div>
--	---------------------------

- Always refer to and follow adhesive manufacturer's complete directions for use for selection and application of the appropriate light cure adhesive. Once the surfaces have been properly treated, they must be kept uncontaminated. Proceed immediately following adhesive application to cementation with Calibra® Veneer Cement, light cured base as indicated in the appropriate section below.

3.4 Veneer cementation – light cured

Before proceeding refer to above sections for Treatment of the Restoration and Tooth Conditioning/Dentin Pretreatment/Adhesive Application.

- Adhesives compatible with light cured, CQ initiated methacrylates are recommended.
- If veneer restorations exceed 1.0mm thickness, are heavily opaqued, or otherwise impede the transmission of light, operator is urged to select Calibra® Esthetic Resin Cement and follow the Dual Cured cementation technique, mixing base with catalyst (available separately).
- Apply adhesive to the internal bonding surface of the restoration if directed by adhesive manufacturer's Directions for Use.

Cementation technique

- Dispense the desired shade of Calibra® Veneer Esthetic Resin Cement base paste from the syringe directly onto the veneer. Protect cement from exposure to light.
- Seat the loaded veneer in place. Remove excess from the gingival margin with a blunt instrument. Briefly (10 seconds or less) light cure the gingival portion only to tack restoration in place. Remove any excess from proximal and lingual margins. **Technique Tip:** Placement of mylar strips between preparation and adjacent teeth prior to seating veneer aids in isolation and in excess cement clean up. After removing gingival excess and tacking veneer in place, remove interproximal excess cement by pulling mylar strip towards the facial, i.e., from tooth to restoration surface.

3.5 Curing

Inadequate polymerization due to insufficient curing

- Check compatibility of curing light.
- Check curing cycle.
- Check curing output before each procedure.

- Using a visible light* curing unit designed to cure CQ initiated methacrylates (spectral output including 470nm), with a minimum output of 550mW/cm², light cure 20 seconds each from the buccal, lingual, and interproximal aspects.
- Following the light cure check and adjust occlusion as necessary. Proceed to Finishing and Polishing section.

*Refer to curing light manufacturer's recommendations for compatibility and curing recommendations.

3.6 Finishing and polishing

- Removal of resin flash is best accomplished with the Enhance® Finishing System of points, cups and discs. The Enhance® System will remove flash and finish restoration margins without removal of or trauma to the enamel.
- Polish final restoration using Enhance® PoGo® Polishing system or Prisma® Gloss® Polishing Paste and Prisma® Gloss® Extra Fine Polishing Paste. (See complete Directions for Use supplied with polishing product chosen)

All shades of Calibra® Veneer Cement are radio-opaque, with a 1mm radiopacity equivalent to 1,77mm radiopacity of aluminum. Aluminum has a radiopacity equivalent to that of dentin. Thus 1mm of material having a radiopacity equivalent to 1mm of aluminum has a radiopacity equivalent to that of dentin.

4. HYGIENE

<div>Cross-contamination</div> <div> <ul style="list-style-type: none">Do not reuse single use products. Dispose in accordance with local regulations. Syringe cannot be reprocessed. Dispose of contaminated syringe in accordance with local regulations.</div>	<div>CAUTION</div>
---	---------------------------

4.1 Cleaning and disinfection

- To prevent syringes from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended to avoid syringe contamination. The use of protective barriers is an additional precautionary measure against gross contamination but not against all contamination.
- Handling with clean/disinfected gloves, syringes may be disinfected with a water-based hospital-level disinfection solution according to national/local regulations.
- Repeated disinfection may damage label. **NOTE:** Destruction of the label by vigorous wiping. Wipe syringe gently.

5. LOT NUMBER AND EXPIRATION DATE

- Do not use after expiration date. ISO standard is used: “YYYY/MM”
- The following numbers should be quoted in all correspondences:
 - Reorder number
 - Lot number
 - Expiration date

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