

CEMENT CLEAN-UP IS A SNAP.

Benefit from the speed of tack-curing without the stress of over-curing. Calibra® Universal self-adhesive resin cement offers a wide tack cure window of up to 10 seconds and an extended 45-second gel phase*, giving you the time you need for a thorough and effective cleanup. It's the perfect cement for nearly any indirect restoration.

Snap up a sample at www.CalibraCement.com and see for yourself why Calibra® Universal is making noise.



Calibra®



It's a Snap.

*10-second tack cure window equals five-second wave cure per surface. For excess cement cleanup, monowave output LED lights with a single peak output around 470nm are recommended. High power, dual or broad spectrum lights may cause premature hardening of excess cement. Check curing light effect on mixed cement in the laboratory prior to clinical use.